

## Curtin MRI Flow Cytometry Facility BD Canto-II

Laser	Fluorescent Channel	LP	Commonly used Fluorescent Dyes
405	450/50 BP		AF405, ATTO 425, BV421, BV480, Pacific Blue, SuperBright 436, Spark Violet 423, StarBright Violet 440, Zombie Violet
	510/50 BP	502	BV510, BV480, Spark Violet 538, StarBright Violet 515, Vio Green
488	530/30 BP	502	FITC, Alexa Fluor™ 488, GFP, CFSE, Fluo-3 ATTO 488, Brilliant Blue 515, CF488A
	585/42 BP	556	PE, PE-CF594, Alexa Fluor 555, PE-eFluor 610, RealBlue 545
	670 LP	655	PC5.5, PC5, PerCP, PerCP-Cy5.5, Brilliant Blue 700, PerCP-eFluor 710, RealBlue 705, PE-Cy7
	780/60 BP	735	PE-Cy7, PE-Fire 810, PE-Vio770, RealBlue 744, RealBlue 780, PE-Alexa 750
640	660/20 BP		APC, Alexa Fluor™ 647, Sytox Red, eFluor™ 660, ATTP 633, Live/Dead Fix Far Red
	780/60 BP	735	APC-A750, APC-Cy7, APC-H7, APC eFluor™ 780, Zombie NIR

<https://app.fluorometer.com/curtin-flowcore>